

Neuro Surgery & Neurology & OrthoPedic

POORVIKA SURGICALS

**#18, 4th Cross, 1st Floor, Bhoopasandra Main Road,
Near Shanimahathma Temple, Nagashettyhalli,
Sanjaynagar Post, Bangalore - 560094
Email : poorvikaot2015@gmail.com**



Niko NeuroSurgery

VYCOR FOR DEEP SEATED BRAIN TUMORS



Niko Expandable Pencil




NIKO SCALP GUN Automatic Preloaded Ranney clip Applicator



Aneurysm Clip Sys
BERSTOCK - GERMANY

NEURO ENDOSCOPIC FORCEPS



MAGIC SUCTION FOR NEUROSURGERY SUCTION AND IRRIGATION



Niko Collagen Glue for Rhinorrhea





Bipolar Forceps & Cables Bayonet / Straight



Spinal Troma Discectomy

AnastoClip AC & GC New –Cleared for RURA Applications Clipped Dural Patch Graft





Raney Clip



NEUROSURGICAL PATTIES


**WHITE
MOP**



Oss Patch



**SIMPST VIDEO
CAMERA CABLE COVER**



DERMIMARKER



Iodrape

External CSF Drainage Systems



**CHHABRA LUMBAR
PERITONEAL SHUNT SYSTEM**



**CHHABRA SLIT N SPRING
HYDROCEPHALUS SHUNT SYSTEM-
VP**



THE VYCOR EXTENSION ARM
FOR USE WITH THE
VYCOR VIEWSITE BRAIN ACCESS SYSTEM
IN ACTION!



Indications :-

- * Adenocarcinoma
- * Brain Tumors



VIEWSITE BRAIN ACCESS SYSTEM (VBAS) 1:1 Ratio

Distal
Part
Size
(with
introducer
removed)



12S

17S

21S

28S


Part No: TC120803

Part No: TC171103

Part No: TC211503

Part No: TC282003

3cm
Working
Length



12M

17M

21M

28M


Part No: TC120805

Part No: TC171105

Part No: TC211505

Part No: TC282005

5cm
Working
Length



12L

17L

21L

28L



Part No: TC120807

Part No: TC171107

Part No: TC211507

Part No: TC282007

7cm
Working
Length



Vycor Medical Inc. 6401 Congress Ave. • Suite 140 • Boca Raton, FL 33487

Phone: 561-558-2020 • Fax: 631-794-2444 • www.vycormedical.com

Aneurysm Clip System

REBSTOCK – GERMANY



STANDARD



TITAN temporär / *TITANIUM temporary*
gold / gold

TITAN permanent / *TITANIUM permanent*
blau / blue

MINI



TITAN temporär / *TITANIUM temporary*
gold / gold

TITAN permanent / *TITANIUM permanent*
pink / pink


STANDARD / MINI



PHYNOX temporär / *PHYNOX temporary*
Federende gold / *Spring gold plated*

PHYNOX permanent / *PHYNOX permanent*
stahlfarben / *steel colored*

- STANDARD-Clips, PHYNOX und *TITAN*
STANDARD-Clips, PHYNOX and *TITANIUM*
- MINI-Clips, PHYNOX und *TITAN*
MINI-Clips, PHYNOX and *TITANIUM*
- Clip-Anlegezangen
Clip Applying Forceps
- Clip-Lagerungstray
Clip Utility Tray




All Sizes both permanent & Temporary
Rebstock – Germany

Anasto Clip for Intra Dural & Extra Dural Procedures

Neurosurgeons

- ❖ In both Intra Dural & Extra Dural Procedures
- ❖ Any procedure involving duraplasty
- ❖ Spinal Trauma
- ❖ Spinal Meningioma (spinal tumor)
- ❖ Craniotomy
- ❖ Posterior Fosse Decompression for Chiari (kee-AR-ee)
- ❖ Malformation
- ❖ Pseudomeningocele (Spine or Skull)
- ❖ Tethered Cord



Orthopedic Spine Surgeons

- ❖ Generally perform Extra Dural Procedures
- ❖ Extra Dural Procedures:
- ❖ Laminectomy / Decompression
- ❖ Discectomy
- ❖ Spine Fusion
- ❖ Hardware removal/revisions
- ❖ Corpectomy (removal of vertebral body)

Anasto Clip Advantages

- ❖ Non Penetrating
- ❖ NO CSF leakage From Suture holes
- ❖ Significantly reduced closure time
- ❖ Ease of use in anatomically restricted areas
- ❖ Immediate hydrostatic
- ❖ Tran nasal Approach



POORVIKA SURGICALS


#18, 4th Cross, 1st Floor, Bhoopasandra Main Road,
Near Shanimahathma Temple, Nagashettyhalli,
Sanjaynagar Post, Bangalore - 560094
Email : poorvikaot2015@gmail.com

Niko- Expandable Cautery with Flexible Tip

(70mm- 200mm)-For Deep Seated Brain Tumors



1. The tip can move to the ideal length, conveniently wise .
2. It can be retracted & extended up to 70 - 200 mm.
3. For All Types of Deep Surgery (Oncology / Onco Surgery / Palvic Surgery / Hip Surgery / CTVS



AnastoClip GC® Closure System

Clipped Dural Patch Graft

No CSF Leakage

Discectomy

Approximation of dural tissue



AnastoClip GC Closure System

Vascular & Neuro

Hardware removal/revisions

Spinal Trauma

Spinal Meningioma

Laminectomy

Craniotomy

Chiari Malformation

Pseudomeningocele

Tethered Cord

Cardial[®]

LeMaitre
VASCULAR

Your Peripheral Vision™

Cardial Surgical Glue

Adhesive Tube Content (g)
15



The cardial surgical glue is a two-component glue consisting of an adhesive component (Gelatin and Resorcinol) and a polymerisation component (Formaldehyde and Glutaraldehyde).

This glue does not contain any blood derivatives and has been developed especially for the use in cardiac and vascular surgery. The glue components are directly mixed at the application site. Providing a specially designed applicator – with Luer-Lock connector – enables an exact application of the polymerisation agent. The polymerisation happens in only 120 seconds.

Typical Applications

- Dissection of vessel layers
- Suture line support with acute aortic dissections
- Support of suture lines during cardiac and vascular procedures

The surgical glue will be delivered ready to use and sterile including the following components in each pack:


- 1 Flexible tube containing the Gelatin Resorcinol mixture
- 1 Opaque glass bottle containing the formaldehyde/glutaraldehyde polymerisation agent
- 1 Luer-Lock single-use syringe, 5 ml
- 1 Applicator with Luer-Lock connector for drawing up and applying the polymerisation agent

POORVIKA SURGICALS

#18, 4th Cross, 1st Floor, Bhoopasandra Main Road,
Near Shanimahatma Temple, Nagashettyhalli,
Sanjaynagar Post, Bangalore - 560094
Email : poorvikaot2015@gmail.com

Minimum Invasive Forceps

Intraoperative Hemostasis Bipolar



Endoscopy Bipolar' —————
Escort for minimally invasive surgery



POORVIKA SURGICALS

#18, 4th Cross, 1st Floor, Bhoopasandra Main Road,
Near Shanimahathma Temple, Nagashettyhalli,
Sanjaynagar Post, Bangalore - 560094
Email : poorvikaot2015@gmail.com

Intraoperative Hemostasis Bipolar Forcep

Application Scope of Endoscopic Neurosurgery

The Neuroendoscopy technology develops a lot in the last 10 years and gradually become mature

Endoscopic Spine Surgery

- *Intraspinal tumor
- Atlanto occipital malformation
- Herniated Disk

ENDOSCOPIC TRANSCRANIAL NEUROSURGERY

- Hemifacial spasm
- Cerebral hemorrhage
- Cholesteatome
- Acoustic neuroma
- Deep Tumor

Endoscopic Ventricular surgery

- Intraventricular tumor ,
- Cerebral cysticercosis
- Hydrocephalus

ENDOSCOPIC TRANSNASAL SKULL BASE SURGERY

- Pituitary Tumor
- MeningiomaAneurysm
- Nasopharyngeal angiofibroma
- Aneurysm

Features

1. Wider View
2. High Intensity Lighting, Strong direct Vision
3. Small Trauma
4. Short Operation Time

PRODUCT ADVANTAGES

1. Elasticity structure to ensure product consistent tensions
2. Ergonomic design handle offering a favorable hand feel
3. Blind Angle
4. Panorama Vision
5. Avoid Bleeding Caused
6. Puncture under

POORVIKA SURGICALS

#18, 4th Cross, 1st Floor, Bhoopasandra Main Road,
Near Shanimahathma Temple, Nagashettyhalli,
Sanjaynagar Post, Bangalore - 560094
Email : poorvikaot2015@gmail.com

NIKO Vascular clamp-Disposable Scalp gun



- Superior to Renny Clips
- No applicator required
- Automatic saves time
- Reusable
- Easy to use
- Very economical & Handy

Model : HJ-II

Functional Features


1. This product is used in neurosurgery to clip the scalp
 2. Convenient for surgical operation
 3. Reduce surgery time
 4. Reduce secondary damage to the human body
 5. Can achieve a seamless clip
- The gun body can hold 15 clips at a time.



Cartridge of 15 Clips

POORVIKA SURGICALS

#18, 4th Cross, 1st Floor, Bhoopasandra Main Road,
Near Shanimahathma Temple, Nagashettyhalli,
Sanjaynagar Post, Bangalore - 560094
Email : poorvikaot2015@gmail.com



NIKO Disposable Skull Driller /Perforator

Functional features

1. Automatic stop is activated when operating on bones with thickness above 3mm.
2. It can be stopped automatically when drilling, providing extra safety.
3. The cutting process is done in one step so that on bone flap will be left over.
4. Circumventions have been made in order to prevent cross infection, risk caused by resterilization and no postoperative cleaning is needed.
5. Reducing the duration of emergency preoperative preparation.
6. Prevention of massive venous sinus bleeding

Design Features

1. Locating cone Preventing skidding
2. Three-blade design: stable and efficient
3. Uniquely designed bezel: extra sharp bezel enhances the speed of drilling
4. One piece, cohesive main and auxiliary drillers: tiny gap between two drillers is designed to protect the automatic stop system from being interrupted by bone fragments.

Model	Specification		
HZ-I	5*8	9*12	11*14
HZ-II	5*8	9*12	11*14
HZ-III	5*8	9*12	11*14
HZ-IV	5*8	9*12	11*14
HZ-V	5*8	9*12	11*14

POORVIKA SURGICALS

#18, 4th Cross, 1st Floor, Bhoopasandra Main Road,
Near Shanimahathma Temple, Nagashettyhalli,
Sanjaynagar Post, Bangalore - 560094
Email : poorvikaot2015@gmail.com

HONGHU
GAOXIANG

Magical Disposable Suction and Washing Head

Exclusive patent counterfeiting inspected

Traditionally, two operators are required when sucking and washing the surface of a wound in surgery. The process can be more difficult for deep operation and endoscopic surgery. By using the Magical Disposable Suction and Washing Head invented by our company, the way of water feeding has been changed. With a very simple switch, the operator can suck or wash the wound at any time. In this way, it can simplify the whole surgical process, improve the quality and reduce the time for a surgery.

Applications: The suction and washing process in neurosurgery, ENT surgery, spinal surgery, endoscopic surgery and other surgeries on deep tissues.



Features:

Stepless regulator design, flushing arbitrary
Special material production, arbitrary bending
Complete specifications, plug and swap easily
Double-sided edge suction, blood spots for suction
Single tube lightweight and delicate, suitable for endoscopic surgery
Sophisticated blunt design, will not damage the tissue
In line with the original habit, left hand right hand operation
Disinfection packaging provided to save patients


Model	External Diameter of the Drainage Catheter	Length
	2±0.5	



1. Automatic stopping with bone thickness above 3 mm
2. It can be stopped automatically when drilling through cranium, providing extra safety
3. No bone flap will be made
4. Prevention of massive venous sinus bleeding.



1. With both single and double retractors, it can be fixed at any point
2. Cross infection is prevented and no cleaning and sterilizing is needed
3. Free angle and optional fixation
4. Highly praised innovation design



1. Very easy to use
2. With quick connected build-in objective lens
3. The objective lens is anti-reflective
4. It is designed with tear tape so that the microscope is protected and there is no need for cutting
5. Infection is avoided

POORVIKA SURGICALS

#18, 4th Cross, 1st Floor, Bhoopasandra Main Road,
Near Shanimahathma Temple, Nagashettyhalli,
Sanjaynagar Post, Bangalore - 560094
Email : poorvikaot2015@gmail.com

NIKO BIPOLAR FORCEPS & CABLES

Reusable / Disposable , European & US Connector ,Straight & Bayonet



Tip Size

- 0.5 mm to 1.5 mm

Length

- 100 mm to 210 mm



SKIN MARKER PEN


New

XL Prep Resistant Ink *colored markers!*

*The gold standard in prep resistant ink skin markers
- now available in colors!*



maximum visibility!



maximum prep resistance!

* results after a 60 second prep with Chlorhex®

Neurology Products

❖ EEG (Electroencephalography)

❖ EMG (Electromyography)

❖ PSG (Polysomnography)

❖ ICU /IOM

A collection of Ten20 Conductive Paste products, including several white jars and boxes, arranged on a light orange background.

Ten20[®] Conductive Paste

Ten20 contains the right balance of adhesiveness and conductivity, enabling the electrodes to remain in place while allowing the transmittance of electrical signals. **Consider the benefits**

Adhesive , Washable and non-drying , Easy clean up , Convenient sizes

A collection of Nuprep Skin Prep Gel products, including several white tubes and boxes, arranged on a light orange background.

4 oz. Tubes (3-pack) , Item #10-30

Nuprep[®] [Nuprep](#) Skin Prep Gel

Nuprep Skin Prep Gel effectively lowers impedance to improve your tracings. Its mild abrasive formula improves conductivity and helps achieve maximum efficiency with your equipment.

A collection of Nuprep Skin Prep Gel products, including a box and several tubes, arranged on a light orange background.

Most effective uses

Audiology
Cardiac catheter monitoring
Cardiac rehabilitation monitoring
ECG telemetry
EEG exams
Evoked potential procedures
Holter monitoring
Sleep tests
Stress tests

A collection of Ten20 Conductive Paste products, including several white jars and boxes, arranged on a light orange background.

2 g Tubes (6-pack)

Item #10-61



Cup Electrodes for EEG, EP and PSG

Cup Electrodes for EEG, EP and PSG

Ref.	Material	Cup Φ	Plug	Cable
0011.014	Sintered silver (Ag/AgCl sint.)	10 mm	1,5 mm female TP	Multicolor L=150cm (60in)
0012.015	Chloride silver (Ag/AgCl)	10 mm	1,5 mm female TP	Multicolor L=150cm (60in)
0012.016	Gold-plated silver (Ag/Au)	10 mm	1,5 mm female TP	Multicolor L=150cm (60in)
0012.017	Silver (Ag)	10 mm	1,5 mm female TP	Multicolor L=150cm (60in)
0012.090	Tin (Sn)	10 mm	1,5 mm female TP	Multicolor L=150cm (60in)



ON REQUEST:

- D=6mm and D=8mm cup for sintered electrodes
- 2mm male Touch Proof plug
- Cable of different lengths (L=100-150-200cm / L=40-60-80in)
- Flat cable from 5 up to 25 electrodes
- Cable of 5 different colours

Key features

- Advanced sensor ensures high signal quality also in long term monitoring (LTM)
- Pre-attached, tangle-free lead wires
- 2 mm hole allows gel to escape during positioning and permits blunt needle insertion
- Molded junction wire-electrode reduces breaking
- Packed individually or in pouches containing 5 units with different colours on cable
- Can be used with most popular gels, pastes and adhesives




Needle Electrodes for EMG and EEG



Concentric Needle Electrodes for EMG

Ref.	Needle Length	Needle Φ	Hub color
0028.012	25mm	0,45mm	Black
0028.013	30mm	0,45mm	Green
0028.014	35mm	0,45mm	Blue
0028.015	40mm	0,45mm	Red
0028.016	45mm	0,45mm	Yellow
0028.018	50mm	0,45mm	Grey
0028.017	65mm	0,55mm	White
0028.112	25mm	0,35mm	Grey
0028.113	30mm	0,35mm	Grey

PREMIUM			
0028.612	25mm	0,45mm	Pink
0028.614	37mm	0,45mm	Orange
0028.618	50mm	0,45mm	Purple



Ask for our
TRIAL KIT !
ref. 0028.K01

Key features

- Stainless steel/precious metal for reliable recordings and low impedance
- Reduced insertion/penetration force for greater patient comfort
- New tip geometry
- Ergonomic colour-coded hub design for ease of use
- Orientation-free connection to BIONEN needle holder
- 100% automated quality control

Needle holder
shielded 5 pole DIN
L=100cm - 40in
EMG system
compatible



Ask for our
TRIAL KIT !
ref. 0028.K02

Ref.	System
0030.012	MYTO 1
0030.013	MICROMED
0030.111	NEMUS NEMUS 2 MYTO 2
0030.112	KEYPOINT

Easin™

Disposable Sclerotherapy Needle

Unique design for easy injection



Easin™

Disposable Sclerotherapy Needle

INDICATIONS

Used to endoscopically introduce a sclerosing agent or vasoconstrictor into selected sites to control actual or potential bleeding lesions in the digestive system; and the injection of saline to aid in Endoscopic Mucosal.



ADVANTAGE FEATURES

- Sterile & Single Use
- Easy Grip Handle
- Handle Clicks on Needle Extension
- Handle Lock Mechanism
- Strong and Smooth Tubing for Frictionless Maneuvering
- Beveled Stainless Steel Needle
- Smooth Metal Tip
- Wide Range of Sizes for Physicians Choice

ORDERING INFORMATION

Type	Needle Length	Needle Inner Diameter	Sheath Length	Outer Diameter
SN18-07/195	4mm	0.7mm (22g)	195cm	2.4mm
SN18-07/230	4mm	0.7mm (22g)	230cm	2.4mm
SN18-11/195	4mm	1.1mm (19g)	195cm	2.4mm
SN18-11/230	4mm	1.1mm (19g)	230cm	2.4mm


POORVIKA SURGICALS

#18, 4th Cross, 1st Floor, Bhoopasandra Main Road,
Near Shanimahathma Temple, Nagashettyhalli,
Sanjaynagar Post, Bangalore - 560094
Email : poorvikaot2015@gmail.com

PREWIRED HEADCAPS

BIONEN prewired headcaps offer:

- Reliable recordings also in long term monitoring (LTM)
- High quality signal
- Molded Electrode with wire chemical corrosion reduced
- Fast connection to the most used EEG systems
- Available customization for montage and electrodes number
- 7 sizes available (1 newborn, 2 pediatric, 4 adult)



SIZES AVAILABLE:

Adult

Cranial Circumference $\pm$ 2cm

XLARGE	60
LARGE	55
MEDIUM	50
SMALL	45

Pediatric


Cranial Circumference $\pm$ 2cm

PEDIATRIC	40
PEDIATRIC SMALL	35
NEWBORN	30



Special headcaps for fMRI and TMS applications and custom montage up to 192 electrodes are available.

Ask for more details!



EAR ELECTRODES

0006.001	Ag/AgCl sintered electrode, cable L=100 cm - 1,5 mm female TP
0006.004	Ag/AgCl sintered electrode, cable L=100 cm - 2 mm male TP
0007.002	Ag/AgCl electrode, cable L=100 cm - 1,5 mm female TP
0007.003	Ag/AgCl electrode, cable L=100 cm - 2 mm male TP
0007.005	Ag/AgCl electrode, cable L=100 cm - 1,5 mm female TP - binauricular
0007.011	Tin electrode, cable L=100 cm - 1,5 mm female TP
0007.012	Tin electrode, cable L=100 cm - 2 mm male TP



POORVIKA SURGICALS

#18, 4th Cross, 1st Floor, Bhoopasandra Main Road,
Near Shanimahatma Temple, Nagashettyhalli,
Sanjaynagar Post, Bangalore - 560094
Email : poorvikaot2015@gmail.com